1. Introduction

The **Ambu® Man** is an advanced instruction and training manikin for the simulation of realistic and correct anatomical and physiological conditions during cardiopulmonary resuscitation by means of rescue breathing and chest compression.

The Ambu Man incorporates a unique, patented hygienic system under which all trainees get their own face piece and head bag, thereby minimising any risk of cross-infection between the trainees. The hygienic system also features realistic expiration from the Ambu Man through the mouth and nose.

With the special hygienic system, disassembly and internal cleaning and disinfection are rendered unnecessary both during and after training.

The Ambu Man is the size of a normal adult and provides an exceptionally lifelike representation of the human anatomy, particularly of those features important to training in modern resuscitation techniques.

The built-in instrumentation shows the effectiveness of the resuscitation, including results achieved for insufflation volume and depth of external chest compression, and indicates any stomach inflation and incorrect hand positioning as well.

For more advanced training and an in-depth analysis of resuscitation effectiveness the Ambu Man can be combined with a computer.

1.1. Differences between models

Your model may not comprise all the features and possibilities mentioned in the directions for use. Please see below for your model.

All models comprise:

1. A hygienic system
2. Carotid pulse
3. Possibility of adjusting chest stiffness

| Model | Monitoring Instrument | Computer Connection | Apex/ sternum Electrodes | Full Body Model |
|--------------------------|-----------------------|---------------------|--------------------------|-----------------|
| B (asic) | | | | /BF |
| C (omputer) | • | • | | /CF |
| D (efibrillation) | • | • | • | /DF |

Models C (torso) and CF (full body) have all the features mentioned except for the possibility of training in defibrillation. The A/D converter with software for computer linkup will have to be purchased as an optional extra.

Models D (torso) and DF (full body) have all the features mentioned above.

2. Specifications

Weight:

Torso with carrying case: approx. 12 kg

Full body with carrying cases: approx. 17 kg

Total length:

Torso: approx. 80 cm

Full body: approx. 170 cm

The guidelines

Ventilation Compression:

AHA/Standard 0,8 - 1,2 l 40 - 50 mm

ERC/Guidelines 2000 with O₂ 0,4 - 0,6 l 40 - 50 mm

Guidelines 2000 without O₂ 0,7 - 1,0 l 40 - 50 mm

Japanese Guidelines 0,5 - 0,8 l 35 - 50 mm

Part/Material

| | |
|---------------------------|------------------------------|
| Basic unit | Polyethylene |
| Instrument part | ABS plastic |
| Skull | PVC, hard |
| Chest skin | PVC, soft |
| Face piece | PVC, soft |
| Head bag, disposable | Polyethylene |
| T-shirt | 100% cotton |
| Tracksuit | 50% cotton and 50% polyester |
| Carrying bag/training mat | PVC coated nylon material |

3. Operating principles - see ①

When air is blown through the mouth or nose of the face piece (1.1), the head bag (1.2) fills, thereby displacing the air in the skull (1.3) through the valve (1.4), which opens if the head is hyperextended or if jaw thrust is being carried out.

The ambient air (not the trainee's air) then passes through the neck tube (1.5) down into the lung bag (1.6) which expands and elevates the chest (1.7) via the spring (1.11).

If insufflation takes place too quickly and/or the air volume is too high, the stomach valve (1.8) will open suddenly and allow air to flow into the stomach bag (1.9), thereby distending the stomach region (1.10). The movements of both lung bag and stomach bag are transmitted mechanically to the monitoring instrument.

When chest compression is performed, the chest and springs (1.11 and 1.12) are squeezed together. The spring (1.12) can be moved to permit changes in chest stiffness. Just as during rescue breathing, the movements during chest compression are transmitted mechanically to the monitoring instrument.

Chest compression must be concluded before rescue breathing is started to reduce the risk of stomach inflation. This applies to Ambu Man as well as to patients.

The carotid pulse can be felt on both sides of the larynx (1.14). During chest compression the pulse is activated automatically by the bulb (1.15), but can be activated manually by the training instructor or by another trainee by squeezing the bellows (1.16).

4. Hygienic system

Schematic diagram:

The head bag has been mounted in the skull and the face piece fitted. Insufflation has just started, see ②.

A: Trainee's air **B:** Ambient air

Insufflation is concluded and expiration has just started, see ③.

NOTE: If a leakage occurs around the mask or if insufflation proves difficult, it is normally due to incorrect positioning when lifting the chin or to insufficient hyperextension of the head.

To comply with the requirements of the hygienic system, and to achieve a proper seal and accurate measuring results, a head bag must always be used.

All trainees must use their own face piece and head bag only during training.

5. Monitoring instrument - see ④

Readings on the AmbuMan monitoring instrument can be viewed from two sides. By activating a cover plate on the side facing the trainee, instrument readings can be concealed from trainees undergoing a CPR test, while allowing the training instructor to monitor the effectiveness of the CPR on the opposite side.

The monitoring instrument is divided into two halves. The left-hand side registers ventilation, including inspiration volume (4.1), and provides indication of stomach inflation (4.2). The right-hand side of the instrument registers chest compression, including depth of compression in millimetres (4.3), and provides indication of incorrect hand positioning (4.4).

The gauges for insufflated volume and depth of compression are designed to display green when the correct volume and chest compression have been reached, and red when the measurement recorded is outside the correct range.

All Ambu training manikins and training products comply with the ILCOR* "Guidelines 2000". Member organizations of ILCOR are AHA, ARC, ERC, HSFC, CLAR, NZRC and RCSA.

Your new Ambu Man has an insufflation volume between 0,7 – 1,0 litres and a depth of compression during chest compression between 40-50 mm for adults.

It is possible to order instrument plates with other recommendations.

The instrument registers inflation in the stomach (4.2) and incorrect hand positioning (4.4) by changing from black to red colour display.

6. Preparation for training

The Ambu Man training manikin is normally supplied in a special carry bag which, when opened, folds out to provide a mat for the trainee during training.

When supplying the Ambu Man as a full body model, the legs will be in a separate carry bag.

6.1. Carry bag/training mat

The strong bag which gives excellent protection to the Ambu Man during transportation can also contain a number of spare face pieces and head bags, see ⑤.

Unfold the bag. Make sure that the manikin is firmly secured to the training mat by inserting the fitting on the training mat into the recess high up on the back of the manikin, see ⑥.

6.2. Mounting the legs on the torso

Unpack the legs.

During transportation the legs can be folded, thus taking up a minimum of space, see ⑦.

Assembly of legs:

Pull down the trousers a little around the hips. Place the legs in a position so that the 2 dowels on the hip part can slide into the corresponding dents at the bottom of the Ambu Man body, see ⑧.

Press the two Velcro straps firmly to the bod, see ⑨.

To remove legs: Pull off the Velcro straps and the legs will be disengaged.

6.3. Mounting the head bag and face piece

Unfold the bag and hold as shown, see ⑩.

Holding the neck of the bag, fold together as shown on the photos, see ⑪ and ⑫.

Insert bag in the head, see ⑬.

Smooth out the edges along the opening in the skull, see ⑭.

Hold the face piece by the ears and position on the skull so that the top is flush with the hairline, see ⑮.

Pull face piece downwards by holding firmly by the ears until in position. Make sure that the edges of the mask are located under the hair at the temples. Fix the openings at the back of the ears on the fixing prongs on each side of the skull, see ⑯.

6.4. Instrument

Activate the instrument by depressing the catch. No damage is done to the Ambu Man nor the instrument during chest compression in case the instrument is not activated (pulled out), see 17.

During the CPR test, the side of the monitoring instrument visible to the trainee may be concealed by the cover plate, see 18.

6.5. Carotid pulse

Fit the tube with bellows onto the tube connector.

The carotid pulse can now be activated manually by the instructor by compressing the bellows, see 19.

6.6. Adjusting chest stiffness

The stiffness of the chest can be adjusted as required by loosening the thumb screw on the back of the manikin: For reduced stiffness set to 'LOW' position, for increased stiffness set to 'HIGH' position, see 20. The values shown, 6N/mm and 11 N/mm (N=Newton), indicate the force which must be applied to compress the chest by 1 mm.

Example: Compressing the chest by 40 mm at the 'LOW' setting, a force of 240 Newton (6Nx40mm) must be applied.

The normal setting is 'MEDIUM' corresponding to approx. 8.5 N/mm.

The Ambu Man can be supplied as a full body model with arms and legs, see 21.

Training may now be started in accordance with the appropriate training programme, see 22 and 23.

NOTE: Remember to change head bags and face pieces for each new trainee.

7. Instructions

The Ambu Man is designed to meet the training recommendations of the International Guidelines 2000 committed to by the American Heart Association (AHA) and the European Resuscitation Council (ERC).

In general, however, it is recommended that instructions should be based on the techniques and provisions applicable in the country in which training will take place.

Ensure there is free access to the airway by hyperextending the head or using jaw thrust manoeuvre, see 24 and 25.

Use the following techniques, see 26:

- Mouth-to-mouth ventilation
- Mouth-to-nose ventilation
- Mouth-to-mask ventilation
- Resuscitator/mask ventilation

The movements of the chest can be clearly observed during insufflation and expiration. The insufflated volume is displayed directly on the monitoring instrument. The expired air can be felt and heard from the mouth and nose of the Ambu Man with no risk of infection, as the air is the same as that insufflated by the trainee into the head bag, see 27.

Accidental stomach inflation is simulated and can be observed in the stomach region and viewed directly on the monitoring instrument.

External chest compression can be practised, and the depth of compression will be displayed on the monitoring instrument in millimetres. It is possible to adjust the stiffness of the chest to simulate a patient with a soft, normal or hard chest. To ensure that the correct point of compression is applied during chest compression, the instrument also provides indication of incorrect hand positioning, see 28.

The carotid pulse can be felt on both sides of the larynx during chest compression, or simulated manually by the instructor, see 29.

8. Removal of face piece and head bag

8.1. Face piece

Pull the ears outwards to release them from the two prongs on each side of the skull. Pull the mask upwards by the ears until free of skull, see 30.

8.2. Head bag

Grip the bag on each side, but avoid closing the opening. Pull bag slowly upwards, twisting from side to side if necessary, to facilitate removal, see 31.

9. Cleaning and disinfection

With the patented hygienic system of the Ambu Man, incorporating exchangeable face piece and head bag, no internal cleaning and disinfection is necessary.

9.1. Head bag

The head bag is disposable and should always be discarded after training.

9.2. Face piece

The face piece may be reused after cleaning and disinfection. Remove the dental insert by pulling it from the face piece as shown in the photos, see 32 and 33.

9.3. Manual cleaning

- A. Rinse face piece and dental insert in clean water.
- B. Wash items in warm water, max. 65 °C (150 °F), using a mild detergent.
- C. Rinse thoroughly in clean water to remove all traces of detergent.

9.4. Machine washing

The face pieces can be washed in an ordinary washing machine. Apply normal dose of detergent and choose a washing programme with maximum temperature of 70 °C (158 °F). In order to avoid the dental inserts rattling against the drum of the washing machine, the face pieces can be placed in a bag made of loosely woven fabric.

9.5. Disinfection

After separating and cleaning the face piece and dental insert, disinfection can be performed as follows:

- A. Place the items in a sodium hypochlorite solution with minimum 500 ppm freely accessible chlorine (1/4 a cup of liquid household bleach per gallon (approx. 4 litres) of tap water for 10 minutes). This solution must be fresh and should be discarded after use.
- B. Place the items in a 70% alcohol chlorhexidine solution for 2 minutes (70% Ethylalcohol and 0.5% chlorhexidine). This method of disinfection is recommended by the Australian Resuscitation Council.
- C. Chemical disinfection can also be performed using recognised disinfectants suitable for use with polyvinylchloride (PVC). The supplier's directions on dosage and disinfection period should be closely followed.
- D. Always rinse the items in clean water after disinfection and allow to dry before storing.
- E. NOTE! The face pieces must not be exposed to boiling, autoclaving or gas sterilisation.

9.6. Cleaning of skull, neck and body

The skull, neck and body of the Ambu Man should be wiped over with a cloth moistened in a mild detergent and then wiped over once more using a cloth moistened in clean water.

When washing, never allow detergent or water to enter the skull, body or instrument unit. If necessary, cover the connection between body and instrument unit with a cloth.

Marks on the Ambu Man left by lipstick or a ball pen can penetrate the material and should therefore be removed as quickly as possible using alcohol.

9.7. Cleaning of clothing

The T-shirt worn by the Ambu Man is made of 100% cotton and is washable at max. 60 °C.

The tracksuit of the full body model is made of 50% cotton and 50% polyester and is washable at max. 40 °C.

9.8. Cleaning of carry bag

The carry bag can be washed in a mild detergent using a cloth or soft brush, rinsed in clean water and then dried.

10. Assembly of face piece

Assemble the face piece by fitting the dental insert as shown and press it into the face piece so that only the flange and collar remain outside, see 32.

First, press the collar of the dental insert into the slot on one side of the face piece.

Next, take the edge and force it over the flange on the dental insert until the flange is seated in the slot in the mask all the way round, see 34.

11. Computer and printer link-up

For more advanced instruction, those models which are equipped for computer link-up can be used together with a computer.

For this purpose Ambu will supply an A/D converter which changes signals obtained from the Ambu Man into signals able to be processed by the computer. Control and calculation software is also supplied for the computer.

Link-up to data-processing equipment enables curve traces to be displayed on the computer monitor during rescue breathing and chest compression.

Further information on the computer program can be found in the Directions for use for the Ambu CPR Software Kit.

Models which are equipped with a computer link-up can also be connected to an Ambu CPR Printer.

Further information about the printer can be found in the Directions for use for the Ambu CPR Printer.

Both the computer and the printer are connected to the Ambu Man via the plug at the side of the instrument.

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