


### Highest Integration of Roles & Skills Showing Real-Time Influence on Patient

This high fidelity manikin provides for life-like interaction of air movement and pressures in the patient. Students feel how chest movement and ventilations can conflict or aid the condition of the patient. It includes state of the art Simulations which show how quality and time critical procedures influence the condition of the patient.

#### Integrated Response for Higher Reality Training



### Smart Debriefing

The After Action Report provides objective data on where performance can be improved. Quick view summaries zoom in on critical areas such as peri-shock analysis and change over of compressors.

### Flexible

Integrates easily into what you are already doing. It is simple to use. Even with two buttons you can get started and receive a valuable report.

If you already have a solution to using your existing AED or Monitor you can continue to use that with SmartMan. If you want to approach it from a new angle, we can provide a total solution.

*Note: Full body manikin, IO Leg, E.T Tube Capable, DTFIS, sensitive to Low Volume, Real Time Feedback Archiving Included, eSERT, Trainer Tools Summary Reports, After Action Report, Simulations, Use Your Existing Solution to Shock, Integrate into Your Existing Megacode Training*

### Simulations

- Real Time Health Status
- Maintenance Threshold
- Responds to Performance
- Shows Optimal Chance for ROSC
- Shows Influence of Fraction

### Tablet Control

- Away from the Direct Action
- Marked in Log and Reports
- Records Critical Events
- Integrates with After Action